



Case Study

Inari Agriculture




Case Study: Inari Agriculture



ESP to supports end to end automation for new seed development and test:

Created workflow management system that handles process flow and data provenance throughout their “High Throughput” Workflow for “biolistics” experiments on their maize and soybean crops. We automate data management in everything from experimental setup, Greenhouse transfers, pollination, harvesting, and reporting; all while providing an intuitive interface for barcode scanning (e.g. scanning plant samples into a 96-well block before sending it off for bioinformatic analysis), inventory management (e.g. for materials and media needed to operate the gene gun), and location management (tracking plant locations on benches in their greenhouses).




Project Scope



This project focused on implementing 8 ESP workflow chains for Inari West.

1. Maize Edit Delivery
2. Soybean Edit Delivery
3. Tomato Edit Delivery
4. Media Lab Request
5. Vectoring
6. Plant Sampling
7. Analysis



Project Scope description

In these workflow chains, users at Inari were able to:

- Pull vector data from Benchling into ESP
- Perform experiment setup and define experiment parameters on which to perform maize edit deliveries
- Enter and track data associated with performing edit delivery experiments
- Generate unique identifiers and/or barcode labels for various entities throughout the experiment, including plants, explants, and treatments
- Perform sampling (by scanning plants into a 96-well plate) and “analysis” (run an automated pipeline to trigger a machine-recommended decision)
- Generate a Greenhouse transfer manifest and track location in the greenhouse
- Establish object hierarchy during pollination

ESP sample ancestry system



Report: Project Sample Tree

Details for Project Sample Tree

Name

Project Sample Tree

Description

Shows sample genealogy for a selected project

Tags

Add a tag

Owner

system admin

Created

Jun 14, 2019 (CDT)

Updated

Jun 14, 2019 (CDT)

(Maximum of 5 generations)

Project

Workflow Testing

Find

(select)

Find

plN000001 (Vector)

EAR000002 (Ear)

MZINW000001 (Plant)

Analysis Result: FAIL / Creation Date: 2019-06-20 15:18:48.928695 / Creator: system admin / Movement Initiator: system admin / Movement Receiver: Inari User / Next Action: Sample & Analyze / Next Step: Proceed to Transfer / Originating Location: Inari-West / plate_name: 96-well Plate 1 / Plate Name: 96-well Plate 1 / Pollinations: (JSON) / Sample Result: FAIL / Self-pollinated?: True / sname: MZINW000001 / Transfer Date: 2019-06-20 15:23:18.192410 / Transfer Location: Inari-West GH1 / Vector: plN000001

BAT000001 (Seed Batch)

Sample ID: BAT000001 / Sample ID Sequence: SEED BATCH SEQUENCE

Samples	Type	Last Modified
▼ TMT000001	Treatment	2019-06-19 10:05:15 AM (CDT)
▼ EAR000002	Ear	2019-06-19 10:06:26 AM (CDT)
▼ MZINW000001	Plant	2019-06-20 10:18:44 AM (CDT)
MZINW000005	Plant	2019-06-20 10:24:08 AM (CDT)
MZINW000006	Plant	2019-06-20 10:24:08 AM (CDT)
MZINW000007	Plant	2019-06-20 10:24:08 AM (CDT)
MZINW000008	Plant	2019-06-20 10:24:08 AM (CDT)
BAT000001	Seed Batch	2019-06-20 10:24:26 AM (CDT)

ESP sample ancestry system keeps track of the relationships between Ears (of corn), explants (a.k.a. embryos harvested from the ears), plants, child plants (as a product of pollination), to batches of seed.

Project management and bottle neck identification: ESP is the master controller of the project



Experiments 9

View and manage all experiments by Workflow

+ New Experiment

+ New Project

Filter Experiments

Advanced Filters [Clear](#)

☐ Hide Completed Experiments

Experiment Name

Status

Project

Workflow

Owner

[By Project](#) [By Workflow](#) [By Owner](#) [List](#)

Workflow : Transformation Setup2 Experiments

Experiment Name	Status	Samples	Project	Owner	Last Activity	
ZM000002	Queued	4	Workflow Testing	system admin	Jun 21, 2019 (CDT)	⋮
ZM000001	Done	4	Workflow Testing	system admin	Jun 19, 2019 (CDT)	⋮

Workflow : Reporting1 Experiments

Experiment Name	Status	Samples	Project	Owner	Last Activity	
ZM000001	Queued	4	Workflow Testing	system admin	Jun 20, 2019 (CDT)	⋮

Workflow : Harvest1 Experiments

Experiment Name	Status	Samples	Project	Owner	Last Activity	
ZM000001	Done	4	Workflow Testing	system admin	Jun 20, 2019 (CDT)	⋮

Workflow : Pollinate1 Experiments

Experiment Name	Status	Samples	Project	Owner	Last Activity	
-----------------	--------	---------	---------	-------	---------------	--

Integration to Benchling ELN to upload vectors created in the lab into ESP



Worksheet: ZM000002
Transformation Setup

Archive

Save

Save and Continue

< Prev

Vector Selection - (4)



Next >

1 of 3 Protocols

[Pull from Benchling](#)






Sample	* Vectors	* Ears	Total Number of Embryos	Notes	<input type="checkbox"/> Complete	≡
<input type="checkbox"/> TMT000004	pIN000001	EAR000004	400		<input type="checkbox"/>	
<input type="checkbox"/> TMT000003	pIN000002	EAR000003	250		<input type="checkbox"/>	
<input type="checkbox"/> TMT000002	pIN000003	EAR000002	370		<input type="checkbox"/>	
<input type="checkbox"/> TMT000001	pIN000004	EAR000001	340		<input type="checkbox"/>	


Columns

Protocol Details

Efficiency report – assess plant viability after treating plants with vector or agrobacterium



LIMS
system admin 

Worksheet: ZM000001
Reporting
Archive
Save

< Prev
Generate Summary Report - (4)
Next >
1 of 1 Protocols


Sample	# Explants	# Plants Harvested	Efficiency (TE%)	Notes	<input type="checkbox"/> Complete	
<input type="checkbox"/> TMT000003	400	250	62%		<input type="checkbox"/>	
<input type="checkbox"/> TMT000002	350	300	85%		<input type="checkbox"/>	
<input type="checkbox"/> TMT000004	250	160	64%		<input type="checkbox"/>	
<input type="checkbox"/> TMT000001	100	20	20%		<input type="checkbox"/>	

Columns
Protocol Details

Process automation: Automatic generation and emailing of plant manifests

ZM000001 | 06/20/2019 | Page 1



Inari Plant Transfer Manifest

Originating Location	Transfer Location
Inari-West	Inari-West GH1
Kurz Purdue Technology Center	3315 Kent Ave.
1281 Win Hentschel Blvd.	West Lafayette, IN 49706
Suite W1108	
West Lafayette, IN 49706	

Plant ID	Date of Transfer	Vector	Creator
MZINW000001	06/20/2019	pIN000001	system admin
MZINW000002	06/20/2019	pIN000001	system admin
MZINW000003	06/20/2019	pIN000001	system admin
MZINW000004	06/20/2019	pIN000001	system admin

Process automation:
Automatic to signal ahead
to the greenhouse what
plants they will be
receiving

Results & Value

L7's ESP enables Inari to meet their goals for data and process management by providing an operations automation + collaboration platform

1. ESP automatically:

- Tracks entity relationships (e.g. Plant → Seed)
- Tracks locations of entities, whether in the lab or on the field
- Validates and processes incoming data

2. Ability to make quick changes to processes in LIMS

3. ESP provides many pieces of automation in the Inari Agriculture implementation. These include:

- Automatic label printing as the user progresses through their workflows
- Calculations to aid in tracking explants
- Integrations with lab equipment (e.g. sequencers, fluorimeters, liquid handlers)

ESP has served Inari Agriculture by managing data collection, accommodating technologies that are currently in use, and accelerating both new and existing processes.